

Star^{ETORQUE} DENTALEZ

POWER, PERFORMANCE AND FLEXIBILITY

Our full-featured electric handpiece system includes a compact and programmable color touchscreen control unit that can be integrated into a variety of delivery units. ETorque offers a broad range of procedural support with multiple attachments that makes it both flexible and efficient — reducing adjustments. Near-instant smooth, quiet power is delivered via the lightweight, brushless motor.

Programmable color touchscreen display

*100-200,000 rpm operating range in
lightweight (62g) brushless motor*

Autoclavable motor and attachments

Stand-alone or delivery unit integration options

Part numbers:

Star ETorque: 265461

Integrated Star ETorque: STAR-ETORQUE



1:5

Friction Grip-L
10,000–200,000 rpm
Cutting, crown & bridge
Part Number: 265471



1:1

AutoLatch-L
1,000–40,000 rpm
Caries removal
Part Number: 265472



5:1

AutoLatch-L
fiber optic light
Part Number: 265475



5:1

AutoLatch
Non-fiber optic
Part Number: 265474



1:1

Straight
1,000–40,000 rpm
Chairside trimming &
denture adjustment
Part Number: 265473



Specifications

Star ETorque & Integrated Star ETorque

Control Unit

Motor Speed Range	Forward/Reverse 100–40,000 rpm
Output Torque	Max 3.5 Ncm
Dimensions	78 mm W x 175 mm H x 100 mm D
System Pressure	3.5–8NL/min
Spray Air	14.5–36.2 psi
Spray Water	11.6–29 psi

Motor

Length	1.25"
Weight	62 grams
dBa	<40



Integrated Star ETorque